

Water Bath

WATER



Compact, ergonomic and ideal for any laboratory bench.

This premium water bath is calibratable and combines precise temperature control with intuitive operation. It is ideal for tempering ointments, emulsions, samples, plates and culture media in laboratory environments, as well as for warm storage and calibration in industrial settings.

PRECISE TEMPERATURE CONTROL AND EXCEPTIONAL FLEXIBILITY

In many laboratory and industrial environments, water baths must deliver uninterrupted reliability even in continuous operation – without taking up too much valuable workspace.

The WTB water bath has been optimised to integrate Memmert's precise and reliable temperature control into a space-saving, robust and durable stainless-steel housing.



KEY PERFORMANCE FEATURES

- 01** Temperature-setting range: +10 to +100 °C
- 02** Calibratable for highly precise temperature control
- 03** Intuitive touchscreen interface for effortless operation
- 04** Compact design for maximising laboratory workspace
- 05** Flexible and expandable with customised accessories



APPLICATIONS

- 01** Tempering of ointments, emulsions, samples, plates and culture media in laboratories
- 02** Warm storage and calibration in industrial settings



PRACTICAL AND ROBUST


Perfect for everyday laboratory use: The high-quality stainless-steel bath with rounded corners is easy to clean and can be completely drained via the standard valve – no tilting required. The housing is manufactured entirely from durable stainless steel too, ensuring an exceptionally long service life. Combined with its compact, ergonomic design, the Memmert WTB is an ideal fit for any laboratory bench.

INTUITIVE AND USER-FRIENDLY

Memmert is renowned for its user-friendly operating concepts – the new WTB water bath is a further testament to this. A particular highlight is the configurable CustomView display, which keeps important data readily accessible at all times.

FLEXIBLE AND MODULAR EXPANSION

Thanks to its modular design, this modern water bath guarantees maximum flexibility. Optional components such as shaking device, insert rack or test-tube holder can be added effortlessly as required.





CONTACT

Learn more about the Memmert WTB and explore detailed technical specifications.

Visit [memmert.com/wtb](https://www.memmert.com/wtb)

For personal consultation, quotations, or technical questions, please contact your local Memmert representative or our international sales team.

Memmert GmbH + Co. KG
P.O. Box 1720 | 91107 Schwabach, Germany
Tel. +49 9122 925-0 | Fax +49 9122 14585
E-Mail: sales@memmert.com
www.memmert.com

 Memmert
 [memmert.family](https://www.facebook.com/memmert.family)
 [memmert.family](https://www.instagram.com/memmert.family)

We reserve the right to make technical changes - dimensions are not binding

Order no. E68754 © by Memmert GmbH + Co. KG